

QUALICHEK™ Aflatoxin M1 Immunoaffinity Column

REF : KBIC1003

Ver 1.0

RUO

Immunoaffinity Column for the use in Quantitative Determination of
Total Aflatoxin M1.

IVD	For Research Use Only	REF	Catalog Number
	Store At	LOT	Batch Code
	Manufactured By		Biological Risk
	Expiry Date		Consult Operating Instructions

For Research Use Only. Purchase does not include or carry the right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of KRISHGEN BioSystems is strictly prohibited.

REF KBIC1003 10 columns

KRISHGEN BioSystems | for US / Europe: toll free +1(888)-970-0827 tel: +1(562)-568-5005
for Asia / India: tel: +91(22)-49198700
email: sales@krishgen.com

Introduction:

Total Aflatoxin M1 is determined in food and feed using Qualichek™ Immunoaffinity Column (Immunoaffinity Chromatography (Qualichek™ IAC)) based on monoclonal antibodies bonded to macroporous pearl cellulose beads in combination with quantification with HPLC (High Performance Liquid Chromatography).

This method of content determination of Aflatoxin M1 s combines the high selectivity of an immunoaffinity column (IAC) with its potential to concentrate eluate with the additional step of purification by HPLC column. If complex matrices are used for testing, many pretreatment methods of Aflatoxin M1 (AFG2, AFG1, AFB2 and AFB1) determination in food and feed will show low sensitivity because of interfering substances.

Intended Use:

The QUALICHEK™ Total Aflatoxin M1 Immunoaffinity Column is used as an analytical tool for determination of Aflatoxin M1 in grains, feed and food preparations.

Specificity

Aflatoxins M1 (AFL) (= AFB1, AFB2, AFG1 and AFG2).

Form:

Columns filled with gel in phosphate buffered saline (50 mM K-phosphate, pH 7.5, 150 mM NaCl) containing 0.05% sodium azide and 0.01% 2-methyl-4-isothiazolin-3-one embedded within two porous frits.

Gel Bed Volume:

0.3 - 0.5ml

Column Working Range:

0.04 - 50 ng AFL.

Lower limit only depends on the sensitivity of the detection system (LC-FLD, LC- MS); the given value represents sensitivity of a common LC-FLD system

Recoveries:

85 - 110% of each Aflatoxins M1 species within the column working range

Organic Solvent Tolerance:

Up to 20% methanol or 10% acetonitrile in loading solution.

pH Tolerance:

6.5 - 8.5

Flow Rate:

1 - 3 ml/min

Elution:

2 - 3 ml Methanol

Column Capacity:

900 ng Aflatoxin M1 s (AFL)

Cross Reactivity:

Aflatoxin B1 (>90%)

Aflatoxin B2 (>90%)

Aflatoxin G1 (>90%)

Aflatoxin G2 (>90%)

Relative recovery rates when a semi-saturative quantity of Aflatoxins M1 (≤ 900 ng) in mixture is applied on the column.

Application:

Clean up of milk samples for reliable identification and accurate quantification of Aflatoxin M1.

Extraction:

- 1) Centrifuge the mlk sample at 1540 g for 15 minutes.
- 2) The cream is then separated using e.g. a separation funnel. If separation is incomplete, re-centrifuge.
- 3) Filter approximately 40-50 ml of sample extract through Whatman No. 4 filter paper immediately.

Note:

The milk sample should be filtered to avoid the fatty particles blocking the immunoaffinity column.

Enrichment Step:

- 1) 50 ml extract is diluted with 5 ml 10X PBS (= 10fold concentrated 50mM PBS buffer, pH=7.4) and then applied in a reservoir on top of the Qualichek™ IAC Aflatoxin M1 column.

To maintain full performance of the column, please make sure that proportion of dilution buffer in the solution on top of the column is not too small.

- 2) Rate of flow through the affinity gel is 1 to 3 ml/min. In case of problematic matrices rate of flow should be below 2ml/min.

Note: Be cautious that there are no big air bubbles in the gel or between gel and the lock outlet of column. This can prevent a permanent flow or the necessary exchange of matter. Depending on application and on expected contents, larger or smaller extract aliquots can be applied. In such cases the sample calculation (see below) must be adapted.

Wash:

- 1) After the entire sample has passed through the gel, wash the gel with 5 ml of 50mM PBS.
- 2) Remove the remaining liquids in the gel by applying either pressure from top of the column or pressure from bottom.

Elution:

- 1) Remove the Sample reservoir on top of the Qualicheck™ IAC Aflatoxin M1 column.
- 2) Place an appropriate vial below the affinity column.
- 3) Elute the bound toxins by using a total of 1.5 ml of methanol as elution solvent. The elution process is performed in two steps to ensure complete release of analytes.
- 4) A volume of 0.5 ml elution solvent is applied. Allow the entire volume to pass through the column.
- 5) Wait for 30 seconds.
- 6) Apply additional 1 ml of elution solvent through the column. The flow rate of the elution process should not exceed 1 ml/min.
- 7) Elute the remaining solvent solutions by application of slight under or over pressure. All methanolic fractions are unified to give the column eluate.

HPLC:

The column eluate may be injected into the HPLC directly.

If the contents are low (<2 ng/gm, related to commodity content), it may be concentrated by evaporation (e.g. using a VLM evaporator at 50°C under a permanent stream of nitrogen).

Note: If the sample is required to be concentrated by evaporation; add 100 ul of acetic acid/water (50/50 v/v) as keeper. To avoid loss of analytes, stop the concentration at a small volume of residue, e.g. 50 to 100 ul. The residue then is redissolved in HPLC solvent (e.g. 0.4ml) and an aliquot is finally injected into the system.

HPLC Conditions:

Instrument: Shimadzu HPLC

Column: Trentec Reprosil Pur RP C18 120 ODS3 5 m; 125x3,0mm with guard column

Mobile Phase A: methanol / deionized water (85/15 v/v)

Mobile Phase B: methanol / acetonitrile / deionized water (18/18/64 v/v/v)

Gradient: 0,01 min B 100 %; 16 min B 100 %; 17 min B 0 %; 19 min B 0 %; 20 min B 100 %

Flow Rate: 0.5ml/min

Time of Analysis: 30min

Injector Volume: 100 ul

Post Column Derivatization

(only advisable in case of concomitant aflatoxin B1 analysis): 32 ppm pyridinium hydrobromide perbromide in dioxan/deionized water (0.1/99.9, v/v)

Flow Rate: 0.2ml/min; reaction capillary put between end of column and detector by T device: PEEK capillary, 1/16'' x 0.25 mm ID; length: 40 cm

Fluorescence Detection: λ EX [nm]: 362nm; λ EM [nm]: 440nm.

Temperature: Machine and eluents are at room temperature. Eluents are degassed with helium gas.

Qualicheck™ IAC Column Characteristics:

A) Working Range and Recovery Rates of Qualicheck™ IAC Aflatoxin M1 Column:

Working Range of Column: 0.04 - 200ng Aflatoxin M1 total per IAC

Zero Contamination of Column: <0.04ng (LOD of HPLC-FLD method)

Guaranteed Recovery Rates (*) within the Working Range:

Aflatoxin M1: >85% +/-5%

* Recovery rates are confined to the complete IAC separation procedure. Recovery rates are confined to solvent content of diluted extract below 20% methanol or 10% acetonitrile.

B) Cross Reactivities of Qualicheck™ IAC Aflatoxin M1 Column:

Aflatoxin M1 (AFM1): 100%

Aflatoxin B1 (AFB1): 98%

Aflatoxin B2 (AFB2): 99%

Aflatoxin G1 (AFG1): 50%

Aflatoxin G2 (AFG2): 98%

** Ratio of IAC recovery rates if a quantity of 5ng Aflatoxin total per column is analyzed

C) Capacity of Qualichek™ IAC Aflatoxin M1 Column:

Maximum Column Capacity: 1.8 ug Aflatoxin M1 total

*** is incubated with the IAC for 5 minutes; then the IAC is washed with 2ml PBS and the non-bonded fraction is analyzed. The difference of added analyte and non-bonded analyte equals maximum column capacity.

Calculation of Results:

Example Sample Calculation of AFM1 content.

A) Calculation of Sample Gram Equivalents per HPLC injection:

$$\frac{50 \text{ gm Sample}}{55 \text{ ml Extraction Solvent}} \times \frac{55 \text{ ml Extract}}{0.4 \text{ ml}} \times \frac{0.1 \text{ ml}}{\text{Injector Volume}} = \frac{12.5 \text{ gm}}{\text{Sample Equivalents}}$$

B) Calculation of Aflatoxin M1 contamination of examined commodity in ng/gm:

$$\frac{\# \text{ ng injected AFM1}}{\text{Sample Equivalents [gm]}} = \text{ng/gm AFM1 in commodity}$$

Storage:

The performance of the immuno affinity column is guaranteed until the expiry date given on the label or certificate if stored at 4°C. Avoid freezing of the product.

Reference:

"Aflatoxin M1 in milk by immunoaffinity column cleanup with TLC/HPLC determination" Luzia Shundo; Myrna Sabino, Brazilian Journal of Microbiology 2006, 37,164-167

LIMITED WARRANTY

Krishgen Biosystems does not warrant against damages or defects arising in shipping or handling, or out of accident or improper or abnormal use of the Products; against defects in products or components not manufactured by Krishgen Biosystems, or against damages resulting from such non-Krishgen Biosystems made products or components. Krishgen Biosystems passes on to customer the warranty it received (if any) from the maker thereof of such non Krishgen made products or components. This warranty also does not apply to Products to which changes or modifications have been made or attempted by persons other than pursuant to written authorization by Krishgen Biosystems.

THIS WARRANTY IS EXCLUSIVE. The sole and exclusive obligation of Krishgen Biosystems shall be to repair or replace the defective Products in the manner and for the period provided above. Krishgen Biosystems shall not have any other obligation with respect to the Products or any part thereof, whether based on contract, tort, and strict liability or otherwise.

Under no circumstances, whether based on this Limited Warranty or otherwise, shall Krishgen Biosystems be liable for incidental, special, or consequential damages.

This Limited Warranty states the entire obligation of Krishgen Biosystems with respect to the Products. If any part of this Limited Warranty is determined to be void or illegal, the remainder shall remain in full force and effect.

Krishgen Biosystems. 2020

THANK YOU FOR USING KRISHGEN PRODUCT!