# SARS-Cov-2 Coronavirus (Covid-19) Real-Time RT-PCR (RT-qPCR) Detection Kit

**REF** : KBRT010

Ver 2.0

**IVD**

SARS-CoV-2 Coronavirus (Covid-19) Real-Time RT-PCR (RT-qPCR) Detection Kit is designed to detect the presence of SARS-CoV-2 Coronavirus in respiratory specimens and serum sample as recommended by CDC (Centre for Disease Control) Atlanta, USA.

<b>IVD</b>	<b>For In-Vitro Diagnostic Use Only</b>	<b>REF</b>	<b>Catalog Number</b>
	<b>Store At</b>	<b>LOT</b>	<b>Batch Code</b>
	<b>Manufactured By</b>		<b>Biological Risk</b>
	<b>Expiry Date</b>		<b>Consult Operating Instructions</b>

*For In-Vitro Diagnostic Use Only. Purchase does not include or carry the right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of KRISHGEN BioSystems is strictly prohibited.*

**KRISHGEN BioSystems**

**REF** KBRT010

 100 reactions

Unit Nos#318/319, Shah & Nahar,  
Off Dr E Moses Road, Worli, Mumbai 400018.  
Tel: 91 (22) 49198700 | Email: sales@krishgen.com

**EC REP**

KinesisDx, Lyoner Strasse 14, Frankfurt, Germany

English Instructions	...	...	...	Pages 3 - 7
French Instructions	...	...	...	Pages 8 - 13
German Instructions	...	...	...	Pages 14 - 19
Spanish Instructions	...	...	...	Pages 20 - 25

**Introduction:**

Coronaviruses are a family of large RNA viruses with size ranging from 26 to 32 kb. These viruses are zoonotic and in human can cause respiratory infections. As the coronavirus is an RNA virus it has a relatively high mutation rate resulting in rapid evolution. In December 2019, a new deadly coronavirus known as SARS-CoV-2 (previously known as 2019-nCoV), which has a high sequence similarity to SARS-CoV, was identified and has caused a pneumonia, known as Covid-19, outbreak in Wuhan, China and spread globally.

**Intended Use:**

SARS-CoV-2 Coronavirus (Covid-19) Real-Time RT-PCR (RT-qPCR) Detection Kit is designed to detect the presence of SARS-CoV-2 Coronavirus in respiratory specimens and serum sample. The kit is designed as per guidelines issued by Centre for Disease Control (CDC), Atlanta, USA.

**Principle:**

This kit contain three primer/probe sets (N gene primer/probe set 1, N gene primer/probe set 2, and N gene primer/probe set 3) that target the conservative regions of coronavirus nucleocapsid (N) gene, the Human RPP30 gene primer/probe set (Human RPP30 gene primer/probe set) that targets exon 1 of human RPP30 gene and serves as a control to assess specimen quality, a non-infectious DNA positive control (Positive control (non-infectious DNA) to ensure reagents and instruments are working properly, and a negative human specimen extraction control (human RNA extract from non-infected samples, Human specimen extraction control) for assessing reverse transcription. Please refer to Tables 4 and 5 for results interpretation.

**Reagents and Materials Provided:**

Sr No	Component	Quantity	Storage
N1	N gene Primer/Probe set 1, lyophilized	1 vial	-20°C
N2	N gene Primer/Probe set 2, lyophilized	1 vial	-20°C
N3	N gene Primer/Probe set 3, lyophilized	1 vial	-20°C
RP	Human RPP30 gene Primer/Probe set, lyophilized	1 vial	-20°C
H2O	Nuclease-free H <sub>2</sub> O	8 ml	4°C
Pos	Positive Control (non-infectious DNA)	200 ul	-80°C
Hsc	Human specimen extraction Control (HSC) (human RNA extract from non-infected samples)	200 ul	-80°C

**Materials to be Provided by the End-User:**

Component	Recommended
RNA samples	user test sample
qPCR plate or tube	
RT-qPCR master mix	TaqPath™ 1-step RT-qPCR Master Mix, CG, 4X (Thermo Fisher)

**Handling / Storage:**

1. Store the Primer/Probes sets, N gene Primer/Probe set 1, N gene Primer/Probe set 2, N gene Primer/Probe set 3 at -20°C in a manual defrost freezer, and the Positive control and Human specimen extraction Control at -80°C, and nuclease-free H<sub>2</sub>O at 4°C.

**Sample Preparation and Storage:**

Wear appropriate personal protective equipment (e.g. gowns, gloves, eye protection) when working with clinical specimens.

Specimen processing should be performed in a certified class II biological safety cabinet following biosafety level 2 or higher guidelines.

1. Respiratory specimens including: nasopharyngeal or oropharyngeal aspirates or washes, nasopharyngeal or oropharyngeal swabs, bronchoalveolar lavage, tracheal aspirates, and sputum.
2. Swab specimens should be collected only on swabs with a synthetic tip (such as polyester or Dacron®) with aluminum or plastic shafts. Swabs with calcium alginate or cotton tips with wooden shafts are not acceptable.

**Assay Procedure:**

1. Prior to use, allow the **Primer/Probe sets (N1, N2, N3 and RP)** to warm to room temperature in the dark.
2. Centrifuge the vials at 1,500 g for 1 minute.
3. Add 200 ul **Nuclease-free H<sub>2</sub>O** to **N gene Primer/Probe set 1 (lyophilized)** and allow to rehydrate for 15 min at room temperature to make N1 primer/probe stock solution. Aliquot as needed. Store at -20°C in a manual defrost freezer. Avoid repeated freeze-and-thaw cycles. Maintain cold and in the dark when thawed.
4. Add 200 ul **Nuclease-free H<sub>2</sub>O** to **N gene Primer/Probe set 2 (lyophilized)** and allow to rehydrate for 15 min at room temperature to make N2 primer/probe stock solution. Aliquot as needed. Store at -20°C in a manual defrost freezer. Avoid repeated freeze-and-thaw cycles. Maintain cold and in the dark when thawed.
5. Add 200 ul **Nuclease-free H<sub>2</sub>O** to **N gene Primer/Probe set 3 (lyophilized)** and allow to rehydrate for 15 min at room temperature to make N3 primer/probe stock solution. Aliquot as needed. Store at -20°C in a manual defrost freezer. Avoid repeated freeze-and-thaw cycles. Maintain cold and in the dark when thawed.
6. Add 200 ul **Nuclease-free H<sub>2</sub>O** to **Human RPP30 gene Primer/Probe set (lyophilized)** and allow to rehydrate for 15 min at room temperature to make RP primer/probe stock solution. Aliquot as needed. Store at -20°C in a manual defrost freezer. Avoid repeated freeze-and-thaw cycles. Maintain cold and in the dark when thawed.
7. With test samples, three control samples should be run concurrently, the **Positive Control (non-infectious DNA)**, the **Human specimen extraction Control (HSC)** and **H<sub>2</sub>O** as the No Template Control (NTC). Prepare four RT-qPCR reactions for each control sample, one with **N1 Primer/Probe stock solution**, one with **N2 Primer/Probe stock solution**, one with **N3 Primer/Probe stock solution**, and one with **RP Primer/Probe stock solution**. Prepare **20 ul RT-qPCR reactions** for one well as shown in Table 1.

**Table 1.**

Control sample	5 ul
Primer/probe stock solution (N1, N2, N3 or	2 ul
1-step RT-qPCR Master mix, 4x	5 ul
Nuclease-free H <sub>2</sub> O	8 ul
<b>Total Volume</b>	<b>20 ul</b>

8. For each extracted RNA test sample, prepare four RT-qPCR reactions, one with **N1 Primer/Probe stock solution**, one with **N2 Primer/Probe stock solution**, one with **N3 Primer/Probe stock solution**, and one with **RP Primer/Probe stock solution**. Prepare **20 ul RT-qPCR reactions** for one well as shown in Table 2.

**Table 2.**

RNA test sample (concentration varies)	5 ul
Primer/Probe stock solution (N1, N2, N3 or RP)	2 ul
1-step RT-qPCR Master mix, 4x	5 ul
Nuclease-free H <sub>2</sub> O	8 ul
<b>Total Volume</b>	<b>20 ul</b>

- Seal the RT-qPCR reaction wells. Centrifuge the plates or tubes at 1,500 g for 15 seconds.
- Setup RT-qPCR reactions as shown in Table 3.


**Table 3.**

Instrument Settings for RT-qPCR reactions. Fluorescence data (FAM) should be collected during the data acquisition step.

Step	Temperature	Time	Number of cycles
UNG Incubation	25°C	2 min	1
Reverse Transcription	50°C	15 min	1
Enzyme Activation	95°C	2 min	1
Denaturation	95°C	3 sec	45
Annealing and Extension	55°C	30 sec	
Data Acquisition	Plate read, detector (FAM)		

**Data Analysis:**

After completion of the PCR run, save and analyze the data following the instrument manufacturer’s instructions. Analysis should be performed separately for each target using a manual threshold setting. Thresholds (Ct Values) should be adjusted to fall within the exponential phase of the fluorescence curves and above any background signal (refer to image below). The procedure chosen for setting the threshold should be used consistently.

**Interpretation of Results:**

**Table 4.** SARS-Cov-2 Kit Control Sample Test Results Interpretation.

Sample	KBRT010-N1	KBRT010-N2	KBRT010-N3	KBRT-RP	Results Interpretation
2020-Pos	+	+	+	+	Expected
	-	-	-	-	PCR Failed
2020-Hsc	-	-	-	+	Expected
	-	-	-	-	Reverse Transcription and/or PCR Failed
	If anyone of three targets is positive			±	Reagent(s) Contaminated
NTC (H <sub>2</sub> O)	-	-	-	-	Expected
	If anyone of four targets is positive				Reagent(s) Contaminated

**Table 5.** CVPD kit target sample test results interpretation when control results are as expected.

KBRT010-N1	KBRT010-N2	KBRT010-N3	KBRT010-RP	Results Interpretation
+	+	+	±	SARS-CoV-2 Detected
If only one, or two, of three targets is positive			±	Inconclusive Result
-	-	-	+	SARS-CoV-2 Not Detected
-	-	-	-	Invalid Result

**Quality Control:**

The primer/probe sets and positive / negative controls are validated by qPCR. The PCR products are analyzed by gel electrophoresis.

**Safety Precautions:**

- **This kit is For In-Vitro Diagnostic Use Only.** Follow the working instructions carefully.
- The expiration dates stated on the kit are to be observed. The same relates to the stability stated for reagents
- Do not use or mix reagents from different lots.
- Do not use reagents from other manufacturers.
- Avoid time shift during pipetting of reagents.
- All reagents should be kept in the original shipping container.
- Some of the reagents contain small amount of sodium azide (< 0.1 % w/w) as preservative. They must not be swallowed or allowed to come into contact with skin or mucosa.
- Source materials maybe derived from human body fluids or organs used in the preparation of this kit were tested and found negative for HBsAg and HIV as well as for HCV antibodies. However, no known test guarantees the absence of such viral agents. Therefore, handle all components and all patient samples as if potentially hazardous.
- Since the kit contains potentially hazardous materials, the following precautions should be observed
  - Do not smoke, eat or drink while handling kit material
  - Always use protective gloves
  - Never pipette material by mouth
  - Wipe up spills promptly, washing the affected surface thoroughly with a decontaminant.
- In any case GLP should be applied with all general and individual regulations to the use of this kit.

**LIMITED WARRANTY**

Krishgen Biosystems does not warrant against damages or defects arising in shipping or handling, or out of accident or improper or abnormal use of the Products; against defects in products or components not manufactured by Krishgen Biosystems, or against damages resulting from such non-Krishgen Biosystems made

products or components. Krishgen Biosystems passes on to customer the warranty it received (if any) from the maker thereof of such non Krishgen made products or components. This warranty also does not apply to Products to which changes or modifications have been made or attempted by persons other than pursuant to written authorization by Krishgen Biosystems.

**THIS WARRANTY IS EXCLUSIVE.** The sole and exclusive obligation of Krishgen Biosystems shall be to repair or replace the defective Products in the manner and for the period provided above. Krishgen Biosystems shall not have any other obligation with respect to the Products or any part thereof, whether based on contract, tort, and strict liability or otherwise. Under no circumstances, whether based on this Limited Warranty or otherwise, shall Krishgen Biosystems be liable for incidental, special, or consequential damages.

This Limited Warranty states the entire obligation of Krishgen Biosystems with respect to the Products. If any part of this Limited Warranty is determined to be void or illegal, the remainder shall remain in full force and effect.

Krishgen Biosystems. 2020

***THANK YOU FOR USING KRISHGEN PRODUCT!***

M/s Genuine Biosystem Private Limited,  
Plot No.97 & 98, Kattabomman Street Parvathi Nagar Extension, Old Perungalathur  
Chennai, Kanchipuram, Tamil Nadu (India) - 600063

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**\*\*\*\*\*DISCLAIMER\*\*\*\*\***

These sequences are intended to be used for the purposes of respiratory virus surveillance and research. The recipient agrees to use them in compliance with all applicable laws and regulations. Every effort has been made to assure the accuracy of the sequences, but CDC cannot provide any warranty regarding their accuracy. The recipient can acknowledge the source of sequences in any oral presentations or written publications concerning the research project by referring to the Division of Viral Diseases, National Center for Immunization and Respiratory Diseases, Centers for Disease Control and Prevention, Atlanta, GA, USA.

**2019-Novel Coronavirus (2019-nCoV) Real-time rRT-PCR Panel****Primers and Probes**

Division of Viral Diseases

<b>2019-Novel Coronavirus (2019-nCoV) Real-time rRT-PCR Panel Primers and Probes</b>				
<b>Name</b>	<b>Description</b>	<b>Oligonucleotide Sequence (5'&gt;3')</b>	<b>Label<sup>1</sup></b>	<b>Working Conc.</b>
2019-nCoV_N1-F	2019-nCoV_N1 Forward Primer	5'-GAC CCC AAA ATC AGC GAA AT-3'	None	20 µM
2019-nCoV_N1-R	2019-nCoV_N1 Reverse Primer	5'-TCT GGT TAC TGC CAG TTG AAT CTG-3'	None	20 µM
2019-nCoV_N1-P	2019-nCoV_N1 Probe	5'-FAM-ACC CCG CAT TAC GTT TGG TGG ACC-BHQ1-3'	FAM, BHQ-1	5 µM
2019-nCoV_N2-F	2019-nCoV_N2 Forward Primer	5'-TTA CAA ACA TTG GCC GCA AA-3'	None	20 µM
2019-nCoV_N2-R	2019-nCoV_N2 Reverse Primer	5'-GCG CGA CAT TCC GAA GAA-3'	None	20 µM
2019-nCoV_N2-P	2019-nCoV_N2 Probe	5'-FAM-ACA ATT TGC CCC CAG CGC TTC AG-BHQ1-3'	FAM, BHQ-1	5 µM
2019-nCoV_N3-F	2019-nCoV_N3 Forward Primer	5'-GGG AGC CTT GAA TAC ACC AAA A-3'	None	20 µM
2019-nCoV_N3-R	2019-nCoV_N3 Reverse Primer	5'-TGT AGC ACG ATT GCA GCA TTG-3'	None	20 µM
2019-nCoV_N3-P	2019-nCoV_N3 Probe	5'-FAM-AYC ACA TTG GCA CCC GCA ATC CTG-BHQ1-3'	FAM, BHQ-1	5 µM
RP-F	RNAse P Forward Primer	5'-AGA TTT GGA CCT GCG AGC G-3'	None	20 µM
RP-R	RNAse P Reverse Primer	5'-GAG CGG CTG TCT CCA CAA GT-3'	None	20 µM
RP-P	RNAse P Probe	5'-FAM – TTC TGA CCT GAA GGC TCT GCG CG – BHQ-1-3'	FAM, BHQ-1	5 µM

<sup>1</sup>TaqMan® probes are labeled at the 5'-end with the reporter molecule 6-carboxyfluorescein (FAM) and with the quencher, Black Hole Quencher 1 (BHQ-1) (Biosearch Technologies, Inc., Novato, CA) at the 3'-end.

Note: Oligonucleotide sequences are subject to future changes as the 2019-Novel Coronavirus evolves.