KRIBIOLISA™ Atezolizumab (Tecentriq™) ELISA

REF : KBI1027

Ver 2.1

IVD

Enzyme Immunoassay for the quantitative determination of Atezolizumab in serum, plasma and cell culture supernatant

IVD	For In-Vitro Diagnostic Use	REF	Catalog Number
	Store At	LOT	Batch Code
	Manufactured By		Biological Risk
	Expiry Date		Consult Operating Instructions

For In-Vitro Diagnostic Use Only. Purchase does not include or carry the right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of KRISHGEN BioSystems is strictly prohibited.

REF KBI1027

 96 tests

KRISHGEN BioSystems | For US / Europe: toll free +1(888)-970-0827 tel: +1(562)-568-5005
For Asia / India: tel: +91(22)-49198700
Email: sales@krishgen.com

EC REP

MDSS GmbH, Schiffgraben 41, 30175 Hannover, Germany

English Instructions	Pages 1 - 13
French Instructions	Pages 14 - 25
German Instructions	Pages 26 - 37
Spanish Instructions	Pages 38 - 49

Introduction:


Atezolizumab (trade name Tecentriq) is a fully humanized, engineered monoclonal antibody of IgG1 isotype against the protein programmed cell death-ligand 1 (PD-L1). In 2015, it was in clinical trials as an immunotherapy for several types of solid tumors. It was under investigation by Genentech/Roche. In April 2016, Roche announced that atezolizumab had been granted fast track status for lung cancer by the FDA. In May 2018, Tecentriq was in combination with Avastin and standard chemotherapy for some patients with lung cancer was granted priority review.

Intended Use:

The KRIBIOLISA™ Atezolizumab ELISA is used as an analytical tool for quantitative determination of Atezolizumab in serum, plasma and cell culture supernatant.

Principle:

The method employs the quantitative sandwich enzyme immunoassay technique. Antibodies to Atezolizumab are pre-coated onto microwells. Samples and standards are pipetted into microwells and human Atezolizumab present in the sample are bound by the capture antibody. Then, a HRP (horseradish peroxidase) conjugated anti-Atezolizumab antibody is pipetted and incubated. After washing microwells in order to remove any non-specific binding, the ready to use substrate solution (TMB) is added to microwells and color develops proportionally to the amount of Atezolizumab in the sample. Color development is then stopped by addition of stop solution. Absorbance is measured at 450 nm.

**Materials Provided:**

1. Anti-Atezolizumab Coated Microtiter Plate (12x8 wells) – 1 no
2. Atezolizumab Standard, (0.5 ml/vial) – 0, 50, 100, 200, 300, 400, 500 and 600 ng/ml
3. Anti-Atezolizumab:HRP Conjugate – 12 ml
4. (2X) Sample Diluent– 25 ml
5. (20X) Wash Buffer– 25 ml
6. TMB Substrate – 12 ml
7. Stop Solution – 12 ml
8. Instruction Manual

Materials to be provided by the End-User:

1. Microtiter Plate Reader able to measure absorbance at 450 nm.
2. Adjustable pipettes and multichannel pipettor to measure volumes ranging from 25ul to 1000ul
3. Deionized (DI) water
4. Wash bottle or automated microplate washer
5. Standard graph paper or software for data analysis
6. Timer
7. Absorbent Paper

Handling/Storage:

1. All reagents should be stored at 2°C to 8°C for stability.
2. All the reagents and wash solutions should be used within 12 months from manufacturing date.
3. Before using, bring all components to room temperature (18-25°C). Upon assay completion ensure all components of the kit are returned to appropriate storage conditions.
4. The Substrate is light-sensitive and should be protected from direct sunlight or UV sources.

Health Hazard Warnings:

1. Reagents that contain preservatives may be harmful if ingested, inhaled or absorbed through the skin.
2. For In-Vitro Diagnostic Use Only.

**Sample Preparation and Storage:**

Blood is taken by venipuncture. Serum is separated after clotting by centrifugation. Plasma can be used, too. Lipaemic, hemolytic or contaminated samples should not be run. Repeated freezing and thawing should be avoided. If samples are to be used for several assays, initially aliquot samples and keep at -20°C.

For Cell Culture Supernatant – If necessary, centrifuge to remove debris prior to analysis. Samples can be stored at -20°C or -80°C. Avoid repeated freeze-thaw cycles.

For Serum / Plasma Test - Samples have to be diluted to 1:1000 (v/v), e.g. 1 ul sample + 999 ul sample diluent prior to assay. The samples may be kept at 2 - 8°C for up to three days. Long-term storage requires -20°C.

Preparation Before Use:

Allow samples to reach room temperature prior to assay. Take care to agitate patient samples gently in order to ensure homogeneity.

Plasma Test Sample preparation - Samples have to be diluted 1:1000 or 1:2000 (v/v), e.g. for 1:1000 (1 ul sample + 999 ul sample diluent) prior to assay. The samples may be kept at 2 - 8°C for up to three days. Long-term storage requires -20°C.

Serum Test Sample preparation - Samples have to be diluted 1:1000 or 1:2000 (v/v), e.g. for 1:2000 (1 ul sample + 1999 ul sample diluent) prior to assay. The samples may be kept at 2 - 8°C for up to three days. Long-term storage requires -20°C.

Reagent Preparation (all reagents should be diluted immediately prior to use):

1. Label any aliquots made with the kit Lot No and Expiration date and store it at appropriate conditions mentioned.
2. Bring all reagents to Room temperature before use.
3. To make Wash Buffer (1X); dilute 50 ml of 20X Wash Buffer in 950 ml of DI water.
4. To make Sample Diluent (1X); dilute 25 ml of 2X Sample Diluent in 25 ml of DI water.

Procedural Notes:

1. In order to achieve good assay reproducibility and sensitivity, proper washing of the plates to remove excess un-reacted reagents is essential.
2. High Dose Hook Effect may be observed in samples with very high concentrations of Atezolizumab. High Dose Hook Effect is due to excess of antibody for very high concentrations of Atezolizumab present in the sample. High Dose Hook effect is most likely encountered from samples early in the purification process. If Hook Effect is possible, the samples to be assayed should be diluted with a compatible diluent. Thus if the Atezolizumab concentration of the undiluted sample is less than the diluted sample, this may be indicative of the Hook Effect.
3. Avoid assay of Samples containing sodium azide (NaN_3), as it could destroy the HRP activity resulting in under-estimation of the amount of Atezolizumab.

4. It is recommended that all Standards and Samples be assayed in duplicates.
5. Maintain a repetitive timing sequence from well to well for all the steps to ensure that the incubation timings are same for each well.
6. If the Substrate has a distinct blue color prior to use it may have been contaminated and use of such substrate can lead to compromisation of the sensitivity of the assay.
7. The plates should be read within 30 minutes after adding the Stop Solution.
8. Make a work list in order to identify the location of Standards and Samples.

Assay Procedure:

1. It is strongly recommended that all Controls and Samples be run in duplicates or triplicates. A standard curve is required for each assay. All steps must be performed at 37°C
2. Add **100 ul** of **Standards** or **Samples** into the respective wells.
3. Cover the plate and incubate for 120 minutes at 37°C
4. Aspirate and wash plate 4 times with **Wash Buffer (1X)** and blot residual buffer by firmly tapping plate upside down on absorbent paper. Wipe of any liquid from the bottom outside of the microtiter wells as any residue can interfere in the reading step.
5. Pipette without delay in the same order **100 ul** of **Anti-Atezolizumab:HRPConjugate** into each well.
6. Cover the plate and incubate for 60 minutes at 37°C
7. Aspirate and wash plate 4 times with **Wash Buffer (1X)** and blot residual buffer by firmly tapping plate upside down on absorbent paper. Wipe of any liquid from the bottom outside of the microtiter wells as any residue can interfere in the reading step.
8. Add **100 ul** of **TMB Substrate** in each well.
9. Incubate the plate at 37°C for 15-30 minutes in dark. DO NOT SHAKE or else it may result in higher backgrounds and worse precision. Positive wells should turn bluish in color.
10. Pipette out **100 ul** of **Stop Solution**. Wells should turn from blue to yellow in color.
11. Read the absorbance at 450 nm with a microplate reader.

Calculation of Results:


Determine the Mean Absorbance for each set of duplicate or triplicate Standards and Samples. Using Standard graph paper, plot the average value (absorbance 450nm) of each standard on the Y-axis versus the corresponding concentration of the standards on the X-axis. Draw the best fit curve through the standard points. To determine the unknown Atezolizumab concentrations, find the unknown's Mean Absorbance value on the Y-axis and draw a horizontal line to the standard curve. At the point of intersection, draw a vertical line to the X-axis and read the Atezolizumab Concentration. If samples were diluted, multiply by the appropriate dilution factor. Software which is able to generate a cubic spline curve-fit is best recommended for automated results.

Note:

It is recommended to repeat the assay at a different dilution factor in the following cases:

- If the sample absorbance value is below the first standard.
- If the absorbance value is equivalent or higher than the 400 ng/ml standard.

Typical Graph

Quality Control:

It is recommended that for each laboratory assay appropriate quality control samples in each run to be used to ensure that all reagents and procedures are correct.

Performance Characteristics of the Kit:

This kit has been validated as per EMA/FDA guidelines in line with ICH Code for Harmonization of Biological Assays.

Sensitivity:

Limit Of Detection: It is defined as the lowest detectable concentration corresponding to a signal of Mean of '0' standard plus 2* SD. 10 replicates of '0' standards were evaluated and the LOD was found to be less than 50 ng/ml.

Specificity:

The antibodies used in the kit are monoclonal antibodies, anti-idiotypic and specific for Atezolizumab. The calibrators / standards used are calibrated against commercially sourced (Tecentriq™).

Linearity:

Standards provided in the kit will be used for measuring the linearity range of Atezolizumab present in matrix.

Precision:

Precision is defined as the percent coefficient of variation (%CV) i.e. standard deviation divided by the mean and multiplied by 100. Assay precision was determined by both intra (n=5 assays) and inter assay (n=5 assays) reproducibility on two pools with low (50ng/ml), medium (300ng/ml) and high (600ng/ml) concentrations. While actual precision may vary from laboratory to laboratory and technician to technician, it is recommended that all operators achieve precision below these design goals before reporting results.

Pool	Intra Assay %CV	Inter Assay %CV
Low	<10%	<10%
Medium	<5%	<5%
High	<5%	<5%

Recovery:

Known amount of Atezolizumab was “spiked” into sample diluent (1X), diluted normal human serum (1:10, 1:100, 1:1000 and 1:2000) and diluted normal human plasma (1:10, 1:100, 1:1000 and 1:2000) and run in the ELISA.

The resulting concentration, or “recovery” of the spiked material, demonstrates if the expected value can be measured accurately. If the recovered value differs significantly from the amount expected, this may be a sign that some factor in the sample matrix may be causing a falsely elevated or falsely depressed value.

It was observed that serum diluted at 1:1000 and 1:2000 and plasma diluted at 1:1000 and 1:2000 obtained the best recoveries (+/-20%).

Atezolizumab (1:10 serum)	Abs1	Abs2	Mean	Interpolated Concentration	Interpolated Concentration % against Actual Concentration
0	0.472	0.481	0.478	86.258	---
50	0.454	0.440	0.444	--	0.00
100	0.481	0.458	0.466	61.714	61.71
200	0.618	0.559	0.579	206.794	103.40
300	0.706	0.703	0.704	310.721	103.57
400	0.793	0.838	0.823	399.475	99.87
500	0.966	0.929	0.941	486.479	97.30
600	1.096	1.097	1.097	607.815	101.30

Atezolizumab (1:100 serum)	Abs1	Abs2	Mean	Interpolated Concentration	Interpolated Concentration % against Actual Concentration
0	0.191	0.194	0.193	---	---
50	0.256	0.236	0.246	64.238	128.48
100	0.270	0.252	0.261	77.821	77.82
200	0.451	0.425	0.438	210.034	105.02
300	0.585	0.570	0.578	300.698	100.23
400	0.762	0.718	0.740	400.205	100.05
500	0.925	0.871	0.898	494.888	98.98
600	1.075	1.084	1.079	602.711	100.45

Atezolizumab (1:1000 serum)	Abs1	Abs2	Mean	Interpolated Concentration	Interpolated Concentration % against Actual Concentration
0	0.095	0.101	0.098	5.766	--
50	0.151	0.140	0.145	55.050	110.10
100	0.185	0.189	0.187	87.435	87.44
200	0.374	0.374	0.374	204.352	102.18
300	0.556	0.564	0.560	305.062	101.69
400	0.744	0.741	0.742	397.973	99.49
500	0.922	0.954	0.938	495.354	99.07
600	1.167	1.145	1.156	602.738	100.46

Atezolizumab (1:2000 serum)	Abs1	Abs2	Mean	Interpolated Concentration	Interpolated Concentration % against Actual Concentration
0	0.072	0.084	0.078	---	---
50	0.107	0.120	0.113	49.984	99.97
100	0.158	0.157	0.157	93.502	93.50
200	0.327	0.315	0.321	211.518	105.76
300	0.510	0.492	0.501	313.836	104.61
400	0.651	0.599	0.625	376.522	94.13
500	0.892	0.915	0.904	503.792	100.76
600	1.144	1.144	1.144	603.071	100.51

Concentration in ng/ml Abs = Absorbance at 450 nm Interpolated Conc = Recovery Concentration post spiking with known standard
 * indeterminate value

Atezolizumab (1:10 plasma)	Abs1	Abs2	Mean	Interpolated Concentration	Interpolated Concentration % against Actual Concentration
0	0.408	0.447	0.434	---	--
50	0.495	0.482	0.486	72.079	144.16
100	0.508	0.499	0.502	92.582	92.58
200	0.582	0.624	0.610	205.576	102.79
300	0.699	0.736	0.724	304.032	101.34
400	0.820	0.866	0.850	401.094	100.27
500	0.923	0.985	0.965	482.913	96.58
600	1.156	1.158	1.157	609.778	101.63

Atezolizumab (1:100 plasma)	Abs1	Abs2	Mean	Interpolated Concentration	Interpolated Concentration % against Actual Concentration
0	0.441	0.467	0.454	0.000	--
50	0.456	0.592	0.524	42.426	84.85
100	0.658	0.647	0.652	116.845	116.85
200	0.837	0.741	0.789	192.876	96.44
300	1.023	0.966	0.994	303.087	101.03
400	1.208	1.169	1.188	404.887	101.22
500	1.315	1.355	1.335	480.861	96.17
600	1.600	1.583	1.591	611.330	101.89

Atezolizumab (1:1000 plasma)	Abs1	Abs2	Mean	Interpolated Concentration	Interpolated Concentration % against Actual Concentration
0	0.115	0.102	0.108	---	---
50	0.126	0.130	0.128	47.925	95.85
100	0.188	0.196	0.192	110.158	110.16
200	0.292	0.397	0.345	196.663	98.33
300	0.562	0.540	0.551	289.076	96.36
400	0.847	0.879	0.863	429.650	107.41
500	0.914	0.971	0.943	469.981	94.00
600	1.196	1.166	1.181	612.714	102.12

Atezolizumab (1:2000 plasma)	Abs1	Abs2	Mean	Interpolated Concentration	Interpolated Concentration % against Actual Concentration
0	0.101	0.101	0.101	---	---
50	0.123	0.120	0.121	44.349	88.70
100	0.190	0.191	0.190	104.677	104.68
200	0.361	0.377	0.369	202.322	101.16
300	0.555	0.593	0.574	293.433	97.81
400	0.835	0.826	0.831	404.221	101.06
500	1.020	1.062	1.041	501.078	100.22
600	1.252	1.206	1.229	598.415	99.74

Concentration in ng/ml Abs = Absorbance at 450 nm Interpolated Conc = Recovery Concentration post spiking with known standard
 * indeterminate value

Limitations of Method

Healthy individuals should be tested negative by the Atezolizumab. Any clinical diagnosis should not be based on the results of in vitro diagnostic methods alone. Physicians are suggested to consider all clinical and laboratory findings possible to state a diagnosis.

Safety Precautions:

- **This kit is for in vitro use only.** Follow the working instructions carefully.
- The expiration dates stated on the kit are to be observed. The same relates to the stability stated for reagents
- Do not use or mix reagents from different lots.
- Do not use reagents from other manufacturers.
- Avoid time shift during pipetting of reagents.
- All reagents should be kept in the original shipping container.
- Some of the reagents contain small amount of sodium azide (< 0.1 % w/w) as preservative. They must not be swallowed or allowed to come into contact with skin or mucosa.
- Source materials maybe derived from human body fluids or organs used in the preparation of this kit were tested and found negative for HBsAg and HIV as well as for HCV antibodies. However, no known test guarantees the absence of such viral agents. Therefore, handle all components and all patient samples as if potentially hazardous.
- Since the kit contains potentially hazardous materials, the following precautions should be observed

- Do not smoke, eat or drink while handling kit material
- Always use protective gloves
- Never pipette material by mouth
- Wipe up spills promptly, washing the affected surface thoroughly with a decontaminant.
- In any case GLP should be applied with all general and individual regulations to the use of this kit.

References:

Phase I dose-finding study of monotherapy with atezolizumab, an engineered immunoglobulin monoclonal antibody targeting PD-L1, in Japanese patients with advanced solid tumors...Hidenori MizugakiNoboru YamamotoEmail authorHaruyasu MurakamiHirotsugu KenmotsuYutaka FujiwaraYoshimasa IshidaTomohisa KawakamiToshiaki Takahashi...Investigational New Drugs...2016...Springer

Bifunctional immune checkpoint-targeted antibody-ligand traps that simultaneously disable TGFβ enhance the efficacy of cancer immunotherapy...Rajani Ravi, Kimberly A. Noonan, Vui Pham, Rishi Bedi, Alex Zhavoronkov, Ivan V. Ozerov, Eugene Makarev, Artem V. Artemov, Piotr T. Wysocki, Raneeh Mehra, Sridhar Nimmagadda, Luigi Marchionni, David Sidransky, Ivan M. Borrello, Evgeny Izumchenko & Atul Bedi...Nature communications...2018...Nature

High expression of PD-L1 in lung cancer may contribute to poor prognosis and tumor cells immune escape through suppressing tumor infiltrating dendritic cells maturation.... Mu CY, Huang JA, Chen Y, Chen C, Zhang XG....Med Oncol....2011...Springer

MPDL3280A (anti-PD-L1) treatment leads to clinical activity in metastatic bladder cancer.... Powles T, Eder JP, Fine GD, Braithen FS, Loriot Y, Cruz C, Bellmunt J, Burris HA, Petrylak DP, Teng SL, Shen X, Boyd Z, Hegde PS, Chen DS, Vogelzang NJ...2014...Nature

Expression of programmed cell death protein 1 (PD-1) and its ligand PD-L1 in colorectal cancer: Relationship with sidedness and prognosis.... Berntsson J, Eberhard J, Nodin B, Leandersson K, Larsson AH, Jirström K... Oncoimmunology...2018....Taylor & Francis


Prognostic impact of programmed cell death-1 (PD-1) and PD-ligand 1 (PD-L1) expression in cancer cells and tumor infiltrating lymphocytes in colorectal cancer....Li Y, Liang L, Dai W, Cai G, Xu Y, Li X, Li Q, Cai S...Mol Cancer....2016...BMC

SCHEMATIC ASSAY PROCEDURE

1. Remove all components, 30 minutes before adding into the assay plate.

2. Avoid repeated cool-thaw of the components as there will be a loss of activity and this can affect the results.

3. Pipette **100 ul Standards / Samples** into the respective wells.

4. Cover plate and incubate for **120 min** at 37°C.

5. Aspirate and wash wells 4 times with **Wash Buffer (1X)**.

6. Pipette **100 ul Anti-Atezolizumab:HRP** into each well.

7. Cover plate and incubate for **60 min** at 37°C.

8. Aspirate and wash wells 4 times with **Wash Buffer (1X)**.

9. Pipette **100 ul TMB Substrate** into each well.

10. Cover plate and incubate for **30 min** at 37°C.

11. Pipette **100 ul Stop Solution** into each well.

12. Read absorbance at 450nm with a microplate reader within **30 min** of stopping reaction.

Typical Example of a Work List

Well #	Contents	Absorbance at 450nm	Mean Absorbance	ng/ml Atezolizumab equivalent
1A 2A	zero std zero std			
1B 2B	50 ng/ml 50 ng/ml			
1C 2C	100 ng/ml 100 ng/ml			
1D 2D	200 ng/ml 200 ng/ml			
1E 2E	300 ng/ml 300 ng/ml			
1F 2F	400 ng/ml 400 ng/ml			
1G 2G	500 ng/ml 500 ng/ml			
1H 2H	600 ng/ml 600 ng/ml			
3A 4A	Sample			
3B 4B	Sample			

LIMITED WARRANTY

Krishgen Biosystems does not warrant against damages or defects arising in shipping or handling, or out of accident or improper or abnormal use of the Products; against defects in products or components not manufactured by Krishgen Biosystems, or against damages resulting from such non-Krishgen Biosystems made products or components. Krishgen Biosystems passes on to customer the warranty it received (if any) from the maker thereof of such non Krishgen made products or components. This warranty also does not apply to Products to which changes or modifications have been made or attempted by persons other than pursuant to written authorization by Krishgen Biosystems.


THIS WARRANTY IS EXCLUSIVE. The sole and exclusive obligation of Krishgen Biosystems shall be to repair or replace the defective Products in the manner and for the period provided above. Krishgen Biosystems shall not have any other obligation with respect to the Products or any part thereof, whether based on contract, tort, and strict liability or otherwise. Under no circumstances, whether based on this Limited Warranty or otherwise, shall Krishgen Biosystems be liable for incidental, special, or consequential damages.

This Limited Warranty states the entire obligation of Krishgen Biosystems with respect to the Products. If any part of this Limited Warranty is determined to be void or illegal, the remainder shall remain in full force and effect.

Krishgen Biosystems. 2020

THANK YOU FOR USING KRISHGEN PRODUCT!


SYMBOLS KEY

	Anti-Atezolizumab Microtiter Plate (12X8 wells)
	Atezolizumab Standard
	Conjugate Horseradish Peroxidase
	(2X) Sample Diluent
	(20X) Wash Buffer
	TMB Substrate
	Stop Solution
	Consult Instructions for Use
	Catalogue Number
	Expiration Date
	Storage Temperature

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