






KRISHZYME™ O-Glycosidase

REF : KPGF-004

Ver 1.0

RIUO

RIUO	For Research & Industrial Use Only	REF	Catalog Number
	Store At	LOT	Batch Code
	Manufactured By		Biological Risk
	Expiry Date		Consult Operating Instructions

For Research and Industrial Use Only. Purchase does not include or carry the right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of KRISHGEN BioSystems is strictly prohibited.

REF KPGF-004


KRISHGEN BioSystems For US/Europe Customers: toll free +1(888)-970-0827 | tel +1(562)-568-5005
For Asia/India Customers: +91(22)-49198700
Email: sales1@krishgen.com | <http://www.krishgen.biz>

Product Description:

KRISHZYME™ O-Glycosidase, also known as Endo- α -N-Acetylgalactosaminidase, specifically catalyzes the removal of Core 1 and Core 3 O-linked disaccharides from glycoproteins.

KRISHZYME™ O-Glycosidase is a recombinant glycosidase cloned from *Enterococcus faecalis* and expressed in *E. coli*.


KRISHZYME™ O-Glycosidase has a molecular weight of 147kDa. This product is a highly specific enzyme, which hydrolyzes the N-acetylgalactosamine glycosidic linkage, liberating the core disaccharide from the serine or threonine residue that is either unsubstituted or present in glycoproteins or glycopeptides.



The enzyme is specific for α -GalNAc linkages, but it has no apparent preference for serine over threonine-linked residues. Any modification of the core structure can block the action of O-glycosidase. If the O-glycan structure is larger than the core structure, for example substituted with N-acetylglucosamine, N-acetylgalactosamine, sialic acid, or fucose, O-glycosidase will not cleave the GalNAc to Ser/Thr linkage.

In such cases the additional monosaccharides must be sequentially hydrolyzed by a series of exoglycosidases, until only the Gal- β (1-3)-GalNAc- or GlcNAc- β (1-3)-GalNAc- core remains. O-Glycosidase can then remove the core structure intact with no modification of the serine or threonine residue.

Usually it is necessary to treat glycoproteins concomitantly with Neuraminidase and O-Glycosidase. Neuraminic Acid residues must be removed in order to allow O-Glycosidase to cleave the O-linked disaccharides. A general Neuraminidase (KPGF-005) or the Protein Deglycosylation-O Kit (for Simple O-Linked Glycans) (KPGF-100) is recommended.



Any modification of the core structure can block the action of O-glycosidase. If the O-glycan structure is larger than the core structure, for example substituted with N-acetylglucosamine, N-acetylgalactosamine, sialic acid, or fucose, O-glycosidase will not cleave the GalNAc to Ser/Thr linkage. In such cases the additional monosaccharides must be sequentially hydrolyzed by a series of exoglycosidases, until only the Gal- β (1-3)-GalNAc- or GlcNAc- β (1-3)-GalNAc-core remains.

Product Size :

Cat No	Pack Size	Concentration
KPGF-004-A	0.07 U / 50 ul	1.4 U /ml
KPGF-004-B	0.28 U / 200 ul	

Physical Form:

KRISHZYME™ O-Glycosidase is supplied as a liquid in 20mM Tris-HCl (pH 7.5 at 25°C), 50mM NaCl and 1mM EDTA at a concentration of 1.4 U/ml.

Reagents Supplied:

The following reagents are supplied with this product:

Composition	Formula	Concentration
Denaturing Buffer	5%SDS, 0.4M DTT	10X
Assay Buffer 2	0.5M Sodium Phosphate (pH7.5 at 25°C)	10X
NP-40 Solution		10%

Product Source:

Recombinant gene Cloned from Enterococcus faecalis and expressed in E.coli.

Product Quality:

≥95% purity, as determined by SDS-PAGE. No other exoglycosidase, endoglycosidase, and protease activity were contaminated.

Unit Definition:

One unit of O-Glycosidase is defined as the amount of enzyme required to remove 0.68 nmol of O-linked disaccharide from 5 mg of neuraminidase digested, non-denatured fetuin in 1 hour at 37°C in a total reaction volume of 100 ul.

Storage Temperature:

Store at -20°C.

Characteristic :

- Glycerol-free for optimal performance in HPLC and mass spectrometry analysis
- ≥95% purity, as determined by SDS-PAGE
- Optimal activity and stability for up to 24months
- Can be used under native or denaturing conditions
- Can be used in conjunction with other exoglycosidases to remove more complex O-glycans
- Recombinant enzyme with no detectable exoglycosidase or other endoglycosidase contaminating activities

Applications:

- Detection of O-glycosylation in proteins
- Epitope and binding site analysis

- Bioactivity in O-glycans
- Analysis of biosynthesis of glycoproteins
- Structural analysis of O-glycan
- Glycoprotein deglycosylation

Suggestions for Use:

Denaturing Reaction Conditions:

- 1) Combine 10-100 ug of glycoprotein, 1 ul of 10X Denaturing Buffer and H₂O (if necessary) to make a 10 ul total reaction volume.
- 2) Denature glycoprotein by heating reaction at 100°C for 10 minutes.
- 3) Make a total reaction volume of 20 ul by adding 2 ul of 10X Assay Buffer 2, 2 ul of 10% NP-40 Solution, 1-2 ul of α 2-3,6,8,9 Neuraminidase, H₂O and 1-2 ul O-Glycosidase, mix gently.
- 4) Incubate reaction at 37°C for 1-4 hours.

Non-Denaturing Reaction Conditions:

- 1) Combine 10-100 ug of glycoprotein, 2 ul of 10X Assay Buffer 2, 1-2 ul of α 2-3,6,8,9 Neuraminidase, H₂O and 1-2 ul O-Glycosidase to make a 20 ul total reaction volume, mix gently.
- 2) Incubate reaction at 37°C for 1-4 hours.

Notes :

- Since O-Glycosidase is inhibited by SDS, it is essential to have NP-40 in the reaction mixture.
- To deglycosylate a native glycoprotein, longer incubation time as well as more enzyme may be required.

References:

Characterization of glycoproteins and their associated oligosaccharides through the use of endoglycosidases. ... F Maley, RB Trimble, AL Tarentino... - Analytical biochemistry, 1989 - Elsevier

Deglycosylation of asparagine-linked glycans by peptide: N-glycosidase F. ... AL Tarentino, CM Gomez, TH Plummer Jr - Biochemistry, 1985 - ACS Publications/

Demonstration of peptide: N-glycosidase F activity in endo-beta-N-acetylglucosaminidase F preparations. ... TH Plummer, JH Elder, S Alexander, AW Phelan... - Journal of Biological ..., 1984 - ASBMB

Glycosylation of *Pichia pastoris* -derived proteins. ... RK Bretthauer, FJ Castellino - Biotechnology and applied ..., 1999 - Wiley Online Library.

Pharmacological chaperones rescue cell-surface expression and function of misfolded V2 vasopressin receptor mutants. ... JP Morello, A Salahpour, A Laperrière... - The Journal of ..., 2000 - Am Soc Clin Investig.

A monolithic PNGase F enzyme microreactor enabling glycan mass mapping of glycoproteins by mass spectrometry. ... AK Palm, MV Novotny - ... An International Journal Devoted to the ..., 2005 - Wiley Online Library.

Protein-protein interaction and not glycosylation determines the binding selectivity of heterodimers between the calcitonin receptor-like receptor and the ... S Hilaiet, SM Foord, FH Marshall, M Bouvier - Journal of Biological ..., 2001 - ASBMB

Purification and Structure-Function Analysis of Native, PNGase F-Treated, and Endo-. beta.-galactosidase-Treated CHIP28 Water Channels. AN Van Hoek, MC Wiener, JM Verbavatz, D Brown... - Biochemistry, 1995 - ACS Publications

Isolation and characterization of CD47 glycoprotein: a multispanning membrane protein which is the same as integrin-associated protein (IAP) and the ovarian tumour ... WJ Mawby, CH Holmes, DJ Anstee, FA Spring... - Biochemical ..., 1994 - biochemj.org

DC-SIGN and L-SIGN are high affinity binding receptors for hepatitis C virus glycoprotein E2 ... S Fong, A Amara, C Houles, F Fieschi... - Journal of Biological ..., 2003 - ASBMB

LIMITED WARRANTY

Krishgen Biosystems does not warrant against damages or defects arising in shipping or handling, or out of accident or improper or abnormal use of the Products; against defects in products or components not manufactured by Krishgen Biosystems, or against damages resulting from such non-Krishgen Biosystems made products or components. Krishgen Biosystems passes on to customer the warranty it received (if any) from the maker thereof of such non Krishgen made products or components. This warranty also does not apply to Products to which changes or modifications have been made or attempted by persons other than pursuant to written authorization by Krishgen Biosystems.

THIS WARRANTY IS EXCLUSIVE. The sole and exclusive obligation of Krishgen Biosystems shall be to repair or replace the defective Products in the manner and for the period provided above. Krishgen Biosystems shall not have any other obligation with respect to the Products or any part thereof, whether based on contract, tort, and strict liability or otherwise.

Under no circumstances, whether based on this Limited Warranty or otherwise, shall Krishgen Biosystems be liable for incidental, special, or consequential damages.

This Limited Warranty states the entire obligation of Krishgen Biosystems with respect to the Products. If any part of this Limited Warranty is determined to be void or illegal, the remainder shall remain in full force and effect.

Krishgen Biosystems. 2025

THANK YOU FOR USING KRISHGEN PRODUCT!

KRISHGEN BIOSYSTEMS®, GENLISA®, DHARMAPLEX™, GENBULK™, GENLISA™, KRISHZYME®, KRISHZYME SIWFT®, KRISHGEN®, KRIBIOLISA®, KRISHPLEX®, TITANIUM®, QUALICHEK® are registered trademarks of KRISHGEN BIOSYSTEMS PRIVATE LIMITED. ©KRISHGEN BIOSYSTEMS. ALL RIGHTS RESERVED.

KRISHGEN BIOSYSTEMS | OUR REAGENTS | YOUR RESEARCH |