

KRISHGEN BIOSYSTEMS

N-Glycan Kit

Simple · Fast · Accurate N-Glycan analysis

**Released, labeled, purified and ready for HILIC-FLR / LC-MS
in under 30 minutes.**



Labeling Chemistry & Total Workflow Time

Labeling method	Label time	Platform	Total analysis
Traditional 2-AB	1–4 h	—	2–3 days
2-AB Express	1 h	Agilent	1.5 h
InstantPC	—	Agilent	1 h
RapiFluor (RFMS)	5 min	Waters / Krishgen	0.5 h



RapiFluor-class speed

A 5-minute rapid-tagging chemistry delivers the fastest release-to-label workflow in the comparison, with full analysis turnaround in roughly half an hour.

Krishgen Fast N-Glycan kit — a cost-effective, performance-equivalent alternative to the Waters RapiFluor kit.

N-linked Glycosylation, decoded

N-glycans are branched sugar chains attached at asparagine (Asn-X-Ser/Thr) residues. On therapeutic antibodies they sit in the Fc CH2 domain and directly shape effector function, half-life and immunogenicity — which is why glycan profiling is a release-critical CQA.

Asn

attachment site


Fc

CH2 glycosylation

CQA

critical quality attr.

Schematic: glycosylated IgG



Fc N-glycan: GlcNAc · Man · Gal · Sia branches

Four Modules, One Streamlined KIT



Release module

PNGase F enzyme + buffer, surfactant for denaturation, and an IgG positive control for process verification.



Labeling module

Rapid fluorescent label reagent with anhydrous DMF for fast, reproducible tagging of released glycans.



Purification module

HILIC-SPE 96-well plate, elution and dilution buffers tuned for high glycan recovery.



Collection module

8-tube strips & lids, collection plate and waste plate for the full release-to-elute flow.

Ordering information

N-Glycan Kit (Single-Column) 24T **KPGF-400-1**

N-Glycan Kit · 24 tests **KPGF-400-2**


N-Glycan Kit · 96 tests **KPGF-400-3**

Sample Prep Done in >30 minutes



Total hands-on prep **< 30 min** before instrument analysis



Comparable Retention Time & Response



What you're seeing


- ✓ Peak retention times overlay within run-to-run tolerance.
- ✓ Relative fluorescence intensities track the reference profile.
- ✓ MS response intensity is comparable across major glycoforms (F(6)A2, F(6)A2G(4)...).

Comparable Retention Time & Response



The retention time of the peak, the fluorescence intensity, and mass spectrum response intensity is comparable to the competitor.

Accurate, well-resolved Glycoform Markers



Marker fidelity



Krishgen IgG positive control profiled on a Waters ACQUITY RDa detector resolves the full G0–G2 glycoform ladder, giving stable retention-time anchors for confident peak assignment across suppliers.

Efficient Release from Native & Recombinant Proteins

Glycoprotein	Source	N-glycan chains
RNase B	Bovine	1
Transferrin	Human	2
Lactoferrin	Human	3
Fetuin	Bovine	3 (+ O-glycans)
IgG	Human	2
GDH (FAD)	P. pastoris	9

Release time-course (10 → 60 min)

SDS-PAGE band shift confirms progressive, near-complete deglycosylation across substrates within one hour.



KEY FEATURES OF OUR KIT

01



Simple & Fast

Complete sample prep in 30 minutes — a cost-effective alternative to the Waters RapiFluor solution.

02



Consistent & Accurate

Comparable fluorescence and mass-spectrometry response intensity versus the reference kit.

03



High Capability

Efficient cleavage of N-glycans from a wide range of native and recombinant glycoproteins.

04



High Flexibility

Scales from 1 to 96 samples — single-column format plus 24- and 96-well plates to choose from.