









KRIBIOLISA™ Leuprolide / Leuprorelin (Lupron) ELISA

REF : KBI5019

Ver.2.1

RUO

Enzyme Immunoassay for the Quantitative Determination of
Leuprolide in human serum and plasma

	For Research Use Only		Catalog Number
	Store At		Batch Code
	Manufactured By		Biological Risk
	Expiry Date		Consult Operating Instructions

For Research Use Only. Purchase does not include or carry the right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of KRISHGEN BioSystems is strictly prohibited.

REF KBI5019

 **96 tests**

Introduction:

Leuprorelin, also known as leuprolide, is a manufactured version of a hormone used to treat prostate cancer, breast cancer, endometriosis, uterine fibroids. Leuprorelin is in the gonadotropin-releasing hormone (GnRH) analogue family of medications. It works by decreasing gonadotropin and therefore decreasing testosterone and estradiol. Leuprorelin was patented in 1973 and approved for medical use in the United States in 1985. It is on the World Health Organization's List of Essential Medicines. It is sold under the brand name Lupron.

Intended Use:

The KRIBIOLISA™ Leuprolide / Leuprorelin (Lupron) is used as an analytical tool for quantitative determination of Leuprolide in human serum and plasma.

Principle:

The Leuprolide ELISA is a competitive immunoassay for the determination of Leuprolide Acetate. Plates are coated with Anti-GnRH antibody. A varying concentration of unlabeled standard or sample and constant concentration of Leuprolide:HRP conjugate will bind in sequence to the antibodies coated on the microplate. Upon washing, unbound Leuprolide : HRP Conjugate will be removed. Bound Leuprolide: HRP Conjugate complex will produce a soluble blue colored product after the addition of TMB Substrate. The enzyme reaction is stopped by dispensing of stop solution into the wells. The optical density (OD) of the solution at 450 nm is inversely proportional to the amount of bound Leuprolide present in the standards or samples.

Materials Provided:

Part	Description	Qty
Anti- GnRH antibody Coated Microtiter Plate	96 well polystyrene microplate (12 strips of 8 wells) coated with Anti-GnRH antibody Coated Microtiter Plate	1 x 96 wells
Leuprolide standard	Leuprolide standard (lyophilized, concentration - 10,000 ng/ml)	2 vials
Leuprolide:HRP Conjugate	Leuprolide:HRP Conjugate prepared in buffer with protein stabilizer and preservatives 0.02% methylisothiazolone and 0.02% bromonitrodioxane.	12 ml
(1X) Sample Diluent	Buffered protein base with preservative thiomersol < 0.01%	50 ml
(1X) Standard Diluent	Buffered protein base with protein stabilizer and preservatives 0.02% methylisothiazolone and 0.02% bromonitrodioxane with 1:10 dilution human serum	10 ml
(20X) Wash Buffer	20-fold concentrated solution of buffered surfactant with preservative thiomersal < 0.01%. May turn yellow over time.	25 ml
TMB Substrate	Stabilized Chromogen	12 ml
Stop Solution	0.73M Phosphoric Acid	12 ml
Instruction Manual		1 no

Materials to be provided by the End-User:

1. Microplate Reader able to measure absorbance at 450 nm.
2. Adjustable pipettes to measure volumes ranging from 50 ul to 1000 ul.
3. Deionized (DI) water.
4. Wash bottle or automated microplate washer.
5. Graph paper or software for data analysis.
6. Tubes to prepare standard/sample dilutions.
7. Timer.
8. Absorbent paper.
9. Incubator

Handling/Storage:

1. All reagents should be stored at 2°C to 8°C for stability.
2. All the reagents and wash solutions should be used within 12 months from manufacturing date.
3. Before using, bring all components to room temperature (18-25°C). Upon assay completion ensure all components of the kit are returned to appropriate storage conditions.
4. The Substrate is light-sensitive and should be protected from direct sunlight or UV sources.

Health Hazard Warnings:

1. Reagents that contain preservatives may be harmful if ingested, inhaled or absorbed through the skin.
2. For Research Use Only.



Sample Preparation and Storage:

Blood is taken by venipuncture. Serum is separated after clotting by centrifugation. Plasma can be used, too. Lipaemic, hemolytic or contaminated samples should not be run. Repeated freezing and thawing should be avoided. If samples are to be used for several assays, initially aliquot samples and keep at -20°C.

For Cell Culture Supernatant – If necessary, centrifuge to remove debris prior to analysis. Samples can be stored at -20°C or -80°C. Avoid repeated freeze-thaw cycles.

For Serum - Samples have to be diluted 1:10 (v/v), e.g. 1 ul sample + 9 ul sample diluent prior to assay. The samples may be kept at 2 - 8°C for up to three days. Long-term storage requires -20°C.

For Plasma - Samples have to be diluted 1:10 (v/v), e.g. 1 ul sample + 9 ul sample diluent prior to assay. The samples may be kept at 2 - 8°C for up to three days. Long-term storage requires -20°C.

Preparation Before Use:

Allow samples to reach room temperature prior to assay. Take care to agitate patient samples gently in order to ensure homogeneity.

Reagent Preparation (all reagents should be diluted immediately prior to use):

1. Label any aliquots made with the kit Lot No and Expiration date and store it at appropriate conditions mentioned.
2. Bring all reagents to Room Temperature before use.
3. To make **(1X) Wash Buffer**; dilute **25 ml of 20X Wash Buffer in 475 ml of DI water**.
4. **Standards Preparation:** Reconstitute the concentrated Standard lyophilized vial with 1 ml of Standard Diluent to obtain 10,000 ng/ml. Keep the vial for 15 mins with gentle agitation and then run the assay procedure. Use the Standard Diluent as the zero standard. Below table shows the calculation for the standard range.

Standard Concentration (ng/ml)	Standard No.	Dilution Particulars
10,000 ng/ml	Lyophilized standard	Lyophilized Standard + 1ml of Standard Diluent
5000 ng/ml	Standard No.6	500ul Reconstituted Standard + 500ul Standard Diluent (1X)
2500 ng/ml	Standard No.5	500ul Standard No.6 + 500ul Standard Diluent (1X)
1250 ng/ml	Standard No.4	500ul Standard No.5 + 500ul Standard Diluent (1X)
625 ng/ml	Standard No.3	500ul Standard No.4 + 500ul Standard Diluent (1X)
312.5 ng/ml	Standard No.2	500ul Standard No.3 + 500ul Standard Diluent (1X)
156 ng/ml	Standard No.1	500ul Standard No.2 + 500ul Standard Diluent (1X)
0	Standard No.0	Only Standard Diluent (1X)

Mix each tube thoroughly before the next transfer. Use the standards for experiment within one hour of preparation of standard. Discard standard after use.

Procedural Notes:

1. In order to achieve good assay reproducibility and sensitivity, proper washing of the plates to remove excess un-reacted reagents is essential.
2. Avoid assay of Samples containing Sodium Azide (NaN₃), as it could destroy the HRP activity of the conjugate resulting in under-estimation of the antibodies.
3. It is recommended that all Standards and Samples be assayed in duplicates or triplicates.
4. Maintain a repetitive timing sequence from well to well for all the steps to ensure that the incubation timings are same for each well.
5. If the Substrate has a distinct blue color prior to use it may have been contaminated and use of such substrate can lead to compromisation of the sensitivity of the assay.
6. The plates should be read within 30 minutes after adding the Stop Solution.
7. Make a work list in order to identify the location of Standards and Samples.
8. Making serial dilution in the wells directly is not permitted.
9. Prepare the Standard within 15 minutes prior to running the assay.
10. Please carefully dilute Standards according to the instruction, and avoid foaming. To minimize imprecision caused by pipetting, use small volumes and ensure that pipettes are calibrated.
11. If crystals have formed in the Wash Solution (20X) concentrate, warm to room temperature and mix gently until the crystals are completely dissolved.
12. Contaminated water or container for reagent preparation will influence the detection results.

Assay Procedure:

1. It is strongly recommended that all Standards and Samples be run in duplicates or triplicates. A standard curve is required for each assay. All steps must be performed at RT.
2. Pipette out **100 ul of Standards and samples** to the respective wells.
3. Add **100 ul Leuprolide:HRP conjugate** to each well.
4. Cover the plate and incubate for 180 mins at RT and shaking at 600rpm.
5. Aspirate and wash plate 4 times with **Wash Buffer (1X)** and blot residual buffer by firmly tapping plate upside down on absorbent paper. Wipe of any liquid from the bottom outside of the microtiter wells as any residue can interfere in the reading step.
6. Add **100 ul of TMB Substrate** in each well.
7. Incubate the plate at RT for 30 minutes in dark. DO NOT SHAKE or else it may result in higher backgrounds and worse precision. Positive wells should turn bluish in color.
8. Pipette out **100 ul of Stop Solution**. Wells should turn from blue to yellow in color.
9. Read the absorbance at 450 nm with a microplate reader.

Calculation of Results:

Determine the Mean Absorbance for each set of duplicate or triplicate Standards and Samples. Using Graph paper, plot the average value (absorbance 450nm) of each standard on the Y-axis versus the corresponding concentration of the standards on the X-axis. Draw the best fit curve through the standard points. To determine the unknown Leuprolide concentrations, find the unknown's Mean Absorbance value on the Y-axis and draw a horizontal line to the standard curve. At the point of intersection, draw a vertical line to the X-axis and read the Leuprolide Concentration. If samples were diluted, multiply by the appropriate dilution factor. Software which is able to generate a cubic spline curve-fit is best recommended for automated results.

Note:


It is recommended to repeat the assay at a different dilution factor in the following cases:

- If the sample absorbance value is below the first standard.
- If the absorbance value is equivalent or higher than the 10000 ng/ml standard.

Typical Data

Standard Concentration (ng/ml)	Mean Absorbance	Interpolated Concentration	% Interpolated Concentration against Actual Concentration
0	2.427	--	--
156	2.389	147.8	94.8
312.5	1.931	351.5	112.5
625	1.698	571.2	91.4
1250	1.297	1442.9	115.4
2500	1.116	2295.1	91.8
5000	0.902	4158.3	83.2
10000	0.571	11775.9	117.8

Typical Graph



abs = absorbance at 450nm

Quality Control:

It is recommended that for each laboratory assay appropriate quality control samples in each run to be used to ensure that all reagents and procedures are correct.

Performance Characteristics:

Please note that this validation is performed in our laboratory and will not necessarily be duplicated in your laboratory. This data has been generated to enable the user to get a preview of the assay and the characteristics of the kit and is generic in nature. We recommend that the user performs at the minimum; the spike and recovery assay and the dilutional linearity assay to assure quality results. For a more comprehensive validation, the user may run the protocols as suggested by us herein below to develop the parameters for quality control to be used with the kit.

Detection Range:

156ng/ml - 10000ng/ml.

Sensitivity:

The minimum detectable dose of Leuprolide is <50 ng/ml.

Precision:

Intra-Assay Precision: 3 samples with low, middle and high level human Leuprolide were tested 20 times on one plate, respectively.

Inter-Assay Precision: 3 samples with low, middle and high level human Leuprolide were tested on 3 different plates, 8 replicates in each plate.

CV (%) = SD/mean X 100

Intra-Assay: CV<15%

Inter-Assay: CV<18%

Safety Precautions:

- **This kit is For Research Use only.** Follow the working instructions carefully.
- The expiration dates stated on the kit are to be observed. The same relates to the stability stated for reagents
- Do not use or mix reagents from different lots.
- Do not use reagents from other manufacturers.
- Avoid time shift during pipetting of reagents.
- All reagents should be kept in the original shipping container.
- Some of the reagents contain small amount of sodium azide (< 0.1 % w/w) as preservative. They must not be swallowed or allowed to come into contact with skin or mucosa.
- Source materials maybe derived from human body fluids or organs used in the preparation of this kit were tested and found negative for HBsAg and HIV as well as for HCV antibodies. However, no known test guarantees the absence of such viral agents. Therefore, handle all components and all patient samples as if potentially hazardous.
- Since the kit contains potentially hazardous materials, the following precautions should be observed
 - Do not smoke, eat or drink while handling kit material
 - Always use protective gloves
 - Never pipette material by mouth
 - Wipe up spills promptly, washing the affected surface thoroughly with a decontaminant.
- In any case GLP should be applied with all general and individual regulations to the use of this kit.




SCHEMATIC ASSAY PROCEDURE

1. Remove all components, 30 minutes before adding into the assay plate.




2. Avoid repeated cool-thaw of the components as there will be a loss of activity and this can affect the results.





3.  Pipette **100 ul Standard and Sample** into respective well.

4.  Pipette **100 ul Leuprolide:HRP conjugate** into the respective wells.

5.  Cover plate and incubate for  **180 mins** RT and shaking at 600 rpm.

6.  Aspirate and wash wells **4 times** with **Wash Buffer (1X)**.

7.  Pipette **100 ul TMB Substrate** into each well.

8.  Cover plate and incubate for  **30 mins** RT.

9. Read absorbance at 450nm with a  microplate reader within  **30 mins** of stopping reaction.

Typical Example of a Work List

Well #	Contents	Absorbance at 450nm	Mean Absorbance	ng/ml Leuprolide
1A	0 ng/ml			
2A	0 ng/ml			
1B	156 ng/ml			
2B	156 ng/ml			
1C	312.5 ng/ml			
2C	312.5 ng/ml			
1D	625 ng/ml			
2D	625 ng/ml			
1E	1250 ng/ml			
2E	1250 ng/ml			
1F	2500 ng/ml			
2F	2500 ng/ml			
1G	5000 ng/ml			
2G	5000 ng/ml			
1H	10000 ng/ml			
2H	10000 ng/ml			
3A	Sample			
4A				
3B	Sample			
4B				

LIMITED WARRANTY

Krishgen Biosystems does not warrant against damages or defects arising in shipping or handling, or out of accident or improper or abnormal use of the product; against defects in products or components not manufactured by Krishgen Biosystems, or against damages resulting from such non-Krishgen Biosystems made products or components. Krishgen Biosystems passes on to customer the warranty it received (if any) from the maker thereof of such non-Krishgen made products or components. This warranty also does not apply to product to which changes or modifications have been made or attempted by persons other than pursuant to written authorization by Krishgen Biosystems.

THIS WARRANTY IS EXCLUSIVE. The sole and exclusive obligation of Krishgen Biosystems shall be to repair or replace the defective product in the manner and for the period provided above. Krishgen Biosystems shall not have any other obligation with respect to the products or any part thereof, whether based on contract, tort, strict liability or otherwise. Under no circumstances, whether based on this Limited Warranty or otherwise, shall Krishgen Biosystems be liable for incidental, special, or consequential damages.

This Limited Warranty states the entire obligation of Krishgen Biosystems with respect to the product. If any part of this Limited Warranty is determined to be void or illegal, the remainder shall remain in full force and effect.













Krishgen Biosystems. 2024

THANK YOU FOR USING KRISHGEN PRODUCT!

KRISHGEN BIOSYSTEMS®, GENLISA®, DHARMAPLEX™, GENBULK™, GENLISA™, KRISHZYME®, KRISHGEN®, KRIBIOLISA®, KRISHPLEX®, TITANIUM®, QUALICHEK® are registered trademarks of KRISHGEN BIOSYSTEMS. ©KRISHGEN BIOSYSTEMS. ALL RIGHTS RESERVED.

KRISHGEN BIOSYSTEMS | OUR REAGENTS | YOUR RESEARCH

SYMBOLS KEY

	Anti-GnRH antibody Coated Microtiter Plate (12x8 wells)
	Standard, lyophilized
	Conjugate Horseradish Peroxidase
	(1X) Sample Diluent
	(1X) Standard Diluent
	(20X) Wash Buffer
	TMB Substrate
	Stop Solution
	Consult Instructions for Use
	Catalog Number
	Expiration Date
	Storage Temperature