

# KRISHZYME™ Recombinant Trypsin

**Catalog Number: KLPL10201 / KLPL10202 / KLP10203**

## Description

Trypsin specifically cleaves C-terminal peptide bonds of lysine and arginine, which can degrade intercellular binding proteins. The amino acid sequence of the recombinant trypsin produced by us is the same as that of porcine pancreas-derived trypsin, and is produced by recombinant Escherichia coli expression. It can replace traditional extraction of trypsin for cell digestion in vaccines, stem cells, immune cell therapy, drug screening, antibodies and other fields. Aprotinin, soybean trypsin inhibitor, etc. can distinctly inhibit enzyme activity. Recombinant trypsin has the same properties as naturally extracted trypsin.

This product is manufactured in an ISO13485 certified laboratory. No animal-origin raw materials are used in the fermentation, purification, and final ingredients.

## Expression Host:

E.coli

## Purity:

>90% as determined by SDS-PAGE quantitative densitometry by Coomassie Blue Staining.  
(RP-HPLC) - Beta-Trypsin >70%  
(RP-HPLC) - Alpha Trypsin < 20%

## Iso-Electric Point:

8.26

## Endotoxin:

< 20 EU/1000 units as determined by the LAL method.

## Host Cell Proteins (E.coli):

< 0.01%

## Host Cell DNA Residue:

< 10ng/ml

## Unit Definition:

25°C, pH 7.6, 3.2ml reaction solution(1cm light path), one trypsin unit (U) was defined as an increase of 0.003 in the absorption value at 253nm by enzymatic hydrolysis of BAEE pe minute.

## Formulation:


KRISHZYME™ Recombinant Trypsin is supplied as a liquid in buffer containing 50 mM Tris-HCl, pH 7.6.

### Molecular Mass:

The KRISHZYME™ Recombinant Trypsin has a calculated molecular mass of ~24 kDa

### SDS-PAGE:

Fig.1.



KDa Marker 10-120

Fig. 1. Purity analysis by SDS-PAGE Detection

### Reconstitution:

Being an enzyme, the concentration may differ from lot to lot. We always recommend referring the accompanying data sheet to view the exact concentration and the recommended dilution schemata.

Centrifuge the vial at 4°C before opening to recover the entire contents. Please contact us for any concerns or special requirements at +91-22-49198700 | Email: sales1@krishgen.com

### Storage:

Store it under sterile conditions at -20°C to -80°C upon receiving for at least 12 months. It is recommended to aliquot the enzyme into smaller quantities for optimal storage. Avoid repeated freeze-thaw cycles.

### Product Instruction

Long-term storage: If the prepared recombinant cell digestion solution needs to be stored for a longer period, Store at -20°C not exceeding 12 months from date of preparation.

### Application:

For Cell Culture

- Tissue block digestion, primary cell acquisition.
- Passage digestion of adherent cells.
- Cell culture by microcarrier method.

- Gently digesting stem cells.
- Immune cell therapy, etc.

For Recombinant Protein

- Recombinant insulin production.
- Protein sequencing, peptide mapping.
- Specific proteolytic processes such as proteomics research.

## References:

Recombinant trypsin production in high cell density fed-batch cultures in Escherichia coli  
L Yee, HW Blanch - Biotechnology and bioengineering, 1993 - Wiley Online Library

Comparison of the effects between animal-derived trypsin and recombinant trypsin on human skin cells proliferation, gene and protein expression  
M Manira, K Khairul Anuar, WT Seet... - Cell and tissue ..., 2014 - Springer

High-level production of recombinant trypsin in transgenic rice cell culture through utilization of an alternative carbon source and recycling system  
NS Kim, HY Yu, ND Chung, TH Kwon... - Enzyme and microbial ..., 2014 - Elsevier

Addressing common challenges of biotherapeutic protein peptide mapping using recombinant trypsin  
T Menneteau, S Saveliev, CI Butré, AKG Rivera... - ... of Pharmaceutical and ..., 2024 - Elsevier

## LIMITED WARRANTY

Krishgen Biosystems does not warrant against damages or defects arising in shipping or handling, or out of accident or improper or abnormal use of the Products; against defects in products or components not manufactured by Krishgen Biosystems, or against damages resulting from such non-Krishgen Biosystems made products or components. Krishgen Biosystems passes on to customer the warranty it received (if any) from the maker thereof of such non Krishgen made products or components. This warranty also does not apply to Products to which changes or modifications have been made or attempted by persons other than pursuant to written authorization by Krishgen Biosystems.

THIS WARRANTY IS EXCLUSIVE. The sole and exclusive obligation of Krishgen Biosystems shall be to repair or replace the defective Products in the manner and for the period provided above. Krishgen Biosystems shall not have any other obligation with respect to the Products or any part thereof, whether based on contract, tort, and strict liability or otherwise. Under no circumstances, whether based on this Limited Warranty or otherwise, shall Krishgen Biosystems be liable for incidental, special, or consequential damages.

This Limited Warranty states the entire obligation of Krishgen Biosystems with respect to the Products. If any part of this Limited Warranty is determined to be void or illegal, the remainder shall remain in full force and effect.

Krishgen Biosystems. 2025

## THANK YOU FOR USING KRISHGEN PRODUCT!

KRISHGEN BIOSYSTEMS®, GENLISA®, DHARMAPLEX™, GENBULK™, GENLISA™, KRISHZYME®, KRISHGEN®, KRIBIOLISA®, KRISHPLEX®, TITANIUM®, QUALICHEK® are registered trademarks of KRISHGEN BIOSYSTEMS. ©KRISHGEN BIOSYSTEMS. ALL RIGHTS RESERVED.

KRISHGEN BIOSYSTEMS | OUR REAGENTS | YOUR RESEARCH |